



SOLIDCLIP™

CA16 Series

10mm
SINGLE USE CLIP APPLIER




Model List

Series / Model No.	CA16 Series							
Description	162	163	164	165	166	167	168	169
Specification	M	M	M	L	L	ML	ML	ML
Pre-load Quantity	15	20	10	10	15	15	20	10
Size of Clip	M: Medium Size			L: Large Size		ML: Medium Large Size		

Clip Size

Drawing / Dimension	Unit: mm		
	W	H	L
M	4.5	4.8	6
ML	4.5	7.8	9
L	4.5	9.8	11



Medscope Biotech Co., Ltd.

Add: 2F, No. 8, Keyi Street, GuangYuan Technology Park,
Zhunan Town, Miaoli County 350 Taiwan, ROC
TEL: +886 37 580055 FAX: +886 37 585050
Website: www.medscope.com.tw